



More than bone regeneration.
A total solution.

More than products.

A comprehensive portfolio.

Specific indications, treatment goals, and patient conditions vary, so our portfolio covers a wide spectrum of solutions, engineered to predictably and reliably repair defects and enhance implant stability.



STRAUMANN® ALLOGRAFT

- Count on this consistently reliable graft material to give your repairs a strong start.
- Delivers strong structural support, rapid bone regeneration, and volume preservation.
- Most similar to a patient's own bone because it's 100% derived from donor bone provided by a world-leading tissue bank.




STRAUMANN® XENOGRAFT*

- With excellent biocompatibility and an optimal balance of calcium and phosphate, this natural substitute has many characteristics comparable to human bone.
- Low crystallinity offers a favorable surface for new bone while high porosity supports osteoconduction.
- A slower resorption rate provides extended volume maintenance.




STRAUMANN® BONECERAMIC™

- Engineered to support restoration and preservation of bone volume, this alloplastic bone substitute is biocompatible and homogenous.
- As it gradually resorbs, it is replaced by new, vital bone.
- Excellent handling characteristics and consistent quality—a dependable performer.



STRAUMANN® MEMBRANE FLEX™*

- Easy to handle, place, and even reposition, this resorbable barrier naturally conforms to defects and contours.
- Its high suture pullout strength means it can be firmly anchored to surrounding tissue.**
- After providing initial stability during wound healing, it predictably resorbs over 3 to 4 months.**



STRAUMANN® MEMBRANE PLUS™*

- This resorbable, crosslinked barrier provides soft tissue support and space maintenance for multiple indications.
- Slow and predictable resorption over 6 months makes it ideal for larger procedures that need more time to heal.
- It's reinforcement that can provide more time for graft stabilization and bone growth.


*May not be available in all countries, subject to regulatory approval and market release. Please ask your Straumann Territory Manager for more information and availability.

**Data on file with manufacturer.

More than options.

Answers to your needs.

Straumann® bone regeneration materials offer predictable remodeling and resorption, with a range of products to meet specific clinical indications and clinician preferences. It's a comprehensive approach to treatment options for oral rehabilitation—all from one provider.

	Indicated	Recommended	Highly Recommended							
										
				Peri-implant/Dehiscence Defect	Fenestration Defect	Sinus Lift	Extraction Socket	Horizontal Augmentation	Furcation Defect	Intraosseous Defect
AlloGraft Ground Cortical			1							
AlloGraft Demineralized Ground Cortical										2
AlloGraft Ground Cancellous										
AlloGraft Cortical/Cancellous Mix										
XenoGraft*										
BoneCeramic™		3				4	5	6		
Membrane Flex™*										
Membrane Plus™*										

Recommendations in this table are based on material characteristics and scientific data (where available). They indicate which Straumann products are best suited to treat the specified defects prior to or after dental implant placement.

*May not be available in all countries, subject to regulatory approval and market release. Please ask your Straumann Territory Manager for more information and availability.


More than a trusted vendor.

A foundation for success.

For 60 years, practitioners like you have trusted Straumann® to deliver products that offer simplicity, precision, and reliability. When you use Straumann as your bone regeneration and implant dentistry vendor, you get even more.

- **High-quality** bone regeneration solutions at competitive prices.
- **Fast, easy, and convenient online ordering**, with 24/7 access to order tracking and time-saving templates for repeat orders.
- **A value-added program** with rewards that grow as your purchases increase—including business development tools, training and education programs, and patient education materials.
- **Best-in-class service and support** that enable you to spend more time with your patients and focus on your priorities.

Work with one knowledgeable Straumann rep—instead of multiple reps from multiple companies—for a streamlined conduit to the products, tools, and answers you need.



More than a dental implant company. A total solution.

When it comes to treatment options, your patients want positive results—both functionally and esthetically. You want reliable solutions that deliver those results.

Straumann®, a leader in dental implant systems, offers a comprehensive and differentiated portfolio for bone regeneration to help satisfy all of these requirements. Plus, with Straumann as your sole source provider, your practice can also benefit from competitive pricing, time-saving convenience, and value-added rewards.

It's a total tooth replacement solution that has the potential to elevate the patient experience and your practice success.



REFERENCES

- ¹ Kolerman R, et al. Histomorphometric analysis of newly formed bone after maxillary sinus floor augmentation using ground cortical bone allograft and internal collagen membrane. J Periodontol. 2008;79(11):2104-2111.
- ² Hoidal MJ, et al. Clinical evaluation of demineralized freeze-dried bone allograft with and without enamel matrix derivative for the treatment of periodontal osseous defects in humans. J Periodontol. 2008;79:2273.
- ³ Schwarz F et al. Guided bone regeneration at dehiscence-type defects using biphasic hydroxyapatite + beta tricalcium phosphate (Bone Ceramic) or a collagen-coated natural bone mineral (Bio-Oss Collagen): an immunohistochemical study in dogs. Int J Oral Maxillofac Surg. 2007;36(12):1198-1206.
- ⁴ Cordaro L et al. Maxillary sinus grafting with Bio-Oss or Straumann Bone Ceramic: histomorphometric results from a randomized controlled multicenter clinical trial. Clin Oral Implants Res. 2008;19(8):796-803.
- ⁵ Mardas N et al. Alveolar ridge preservation with guided bone regeneration and a synthetic bone substitute or a bovine-derived xenograft: a randomized, controlled clinical trial. Clin Oral Impl Res. 2010;21(688-698).
- ⁶ Mardas N et al. Radiographic alveolar bone changes following ridge preservation with two different biomaterials. Clin Oral Implants Res. 2011;22(4):416-423.

To learn more, **contact your Straumann Territory Manager**
or contact Customer Service.

Buy online at straumann.us/eShop.



International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone: +41 (0)61 965 11 11
Fax: +41 (0)61 965 11 01

Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone: 800.448.8168 (US) • 800.363.4024 (CA) Fax:
978.747.2490
www.straumann.us • www.straumann.ca