

Guiding your Growth
with Straumann® Guided Surgery
Flexible Guided Surgery Kit Solutions

Four Straumann® Guided Surgery Kits


- Instrument Kits designed to give you a choice of workflow options
- Distinct color coded diameter grommets
- Auxiliary instruments organized in the lower tray
- Fully sterilizable cassettes

Comprehensive BL/TL Guided Surgery Kit

The Comprehensive BL/TL Guided Kit contains the following:

- Complete set of BL and TL Guided Surgery Drills to prepare any osteotomy
- Milling Cutters for preparing the alveolar ridge
- Guided Surgery Drill Handles
- Taps and Profile Drills and C Handles
- Template Fixation Pins
- Complete set of drivers and torque wrench

NOTE: Fully Guided Surgery for BL and TL Implants requires our Comprehensive Guided Surgery Kit. A coDiagnostiX® Treatment Plan is required.



Comprehensive BL/TL Guided Kit

The Straumann® Comprehensive Guided Surgery Kit Supports:

- Experienced Surgeons
- The Synergy Workflow
 - Fully Guided Surgery including Guided Implant placement
 - Pre-Operative Provisionalization
 - Tissue Management
- Bone, Tissue or Tooth Supported Protocols
- Flapless Surgery

NOTE: Owning coDiagnostiX® treatment planning software is not required. Many Straumann Alliance Labs or Service Provides are available to assist you with your treatment planning requirements.

NOTE: Bone Level Tapered Instruments can be purchased separately: USKGSBLTREPL in US; CAKGSBLTREPL in Canada

| | Comprehensive BL/TL Kit Order Number |
|--------|---|
| US | USKGS-C |
| Canada | CASET CGSK |

Basic BL or BLT Guided Surgery Kits

For Bone Level or Bone Level Tapered Implants, each respective kit contains the following:

- Complete set of Guided Surgery Drills to prepare any osteotomy for Bone Level or Bone Level Tapered Implants
- Milling Cutters for preparing the alveolar ridge
- Guided Surgery Drill Handles
- Taps and Profile Drills and C Handles
- Template Fixation Pins
- Complete set of drivers and torque wrench

NOTE: The BL or BLT Basic Guided Surgery Kit is specific for Bone Level or Bone Level Tapered Implants. Fully Guided Surgery can be performed. A coDiagnostiX® Treatment Plan is required.



BL or BLT Basic Guided Kit

The Basic Guided Surgery Kits Support:

- Experienced Surgeons placing BL or BLT only
- The Synergy Workflow
 - Fully Guided Surgery including Guided Implants
 - Pre-Operative Provisionalization
 - Soft Tissue Management
- Bone, Tissue or Tooth Supported Protocols
- Flapless Surgery
- Drilling protocols from coDiagnostiX® Treatment Plans

NOTE: Owning coDiagnostiX® treatment planning software is not required. Many Straumann Alliance Labs or Service Providers are available to assist you with your treatment planning requirements.

| | BL Kit Order Number | BLT Kit Order Number |
|--------|------------------------|-------------------------|
| us | USKGSBLBASIC | USKGSBLTBASIC |
| Canada | CAKGSBLBASIC | CAKGSBLTBASIC |

Guided Bed Preparation for BL and BLT Kits

A Flexible Solution for Basic Guided Bed Preparation Only


A cost effective solution for Straumann® Surgeons
Available for Bone Level or Bone Level Tapered Implants

- Complete set of Guided Surgery Drills to prepare any osteotomy
- Guided Surgery Drill Handles
- Guided Taps or Profile Drills Not Included

The Bed Preparation Guided Surgery Kits Support:

- Straumann Surgeons
- Implant Bed Preparation Only
- Drilling protocols from coDiagnostiX® Treatment Plans

IMPORTANT: A Guided Bed Preparation Kit allows the Surgeon to prepare the osteotomy and the implant bed. The remaining steps for tapping and profiling are carried out with conventional Straumann Instruments.



| | BL Kit Order Number | BLT Kit Order Number |
|--------|------------------------|-------------------------|
| us | USKGS-BLBASIC | USKGSBLTBBED |
| Canada | CAKGSBBED | CAKGSBLTBBED |

Straumann® Guided Pilot Surgery Kits

A Flexible Solution for Pilot Surgical Guides

Starting with the **2.2 mm Pilot Drill** through a 2.2 mm Guided Sleeve

Contains:

- Three 2.2 mm Pilot Drills in three lengths
- One Mini Surgical Cassette



Basic Guided Pilot Kit - 2.8 mm

Starting with the **2.8 mm Twist Drill Pro** through a 2.8 mm Guided Sleeve

Contains:

- Three 2.8 mm Twist Drill Pro in three lengths
- One Milling Cutter 2.8 mm
- One Fixation Pin
- One Mini Surgical Cassette

The Straumann® Guided Pilot Surgery Kits Support:

- Straumann Surgeons
- Surgeons New to Guided Surgery
- Pilot drilling protocols from coDiagnostiX® Treatment Plans

IMPORTANT: A Pilot Surgical Guide is used to drill the first step of the osteotomy. A coDiagnostiX® Treatment Plan is required. The remaining osteotomy steps are carried out with conventional Straumann Instruments.

NOTE: Owning coDiagnostiX® treatment planning software is not required. Many Straumann® Alliance Labs or Service Provides are available to assist you with your treatment planning requirements.

| | 2.2 mm BL/TL Kit Order Number | 2.2 mm BLT Kit Order Number |
|--------|----------------------------------|--------------------------------|
| US | USKBLPILOT2 | USKBLTPILOT2 |
| Canada | CAKBLPILOT2 | CAKBLTPILOT2 |

| | 2.8 mm BL/TL Kit Order Number | 2.8 mm BLT Kit Order Number |
|--------|----------------------------------|--------------------------------|
| US | USKBLPILOT1 | USKBLTPILOT1 |
| Canada | CAKBLPILOT1 | CAKBLTPILOT1 |

Don't see the kit you're looking for?

We're committed to providing dental professionals with the proper kits to facilitate guided surgery. Contact your local Straumann Territory Manager for assistance with kit customization.

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01

Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168 (US) • 800/363 4024 (CA)
Fax 978/747 2490
www.straumann.us • www.straumann.ca

© Straumann USA, LLC 2016. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.