



SIMPLY BEAUTIFUL RESTORATIONS



COMMITTED TO
SIMPLY DOING MORE
FOR DENTAL PROFESSIONALS

RESTORE **SIMPLY**



Simplified and secure connections

- CrossFit™ Connection to Straumann® bone level implants provides:
 - Stability through deep conical connection and patented grooves
 - Clear path of insertion
- Morse taper connection to Straumann tissue level implants provides a secure frictional fit



Solid Abutment all-in-one sets

- Helps make restoring easy by packaging all compatible components – solid abutment, positioning cylinder, impression cap, analog, plastic coping, temporary coping – in a single kit

One restorative kit for various indications

- Whether you're doing – single-crown, multiple unit, or overdentures – there's one restorative kit

Prosthetic planning kits

- Plan the ideal esthetics for your patients with prosthetic planning kits for tissue level and bone level implants



Platform flexibility

- A choice of abutment heights, diameters, angulation, and materials are designed to create the optimal fit for virtually any case, creating a beautiful esthetic result


Bone Level						
	Cementable Prosthetic handling closest to a natural tooth preparation	Anatomic Esthetic cement-retained crowns and bridges	Multi-Base Flexible base for screw-retained multiple unit restorations	Gold Versatile solution for a wide range of uses	Customized Straumann customized implant prosthetics	LOCATOR® bar Hybrid solutions for edentulous cases
Tissue Level						
	Solid Prosthetic handling closest to a natural tooth preparation	synOcta® line Flexible crown and bridge restorations		Customized Straumann customized implant prosthetics	LOCATOR, bar, retentive anchor Hybrid solutions for edentulous cases	

LOCATOR® is a registered trademark of Zest Anchors, Inc.

RESTORE **BEAUTIFULLY**

Simplified soft tissue management

- Whether you restore tissue level or bone level implants, Straumann® tissue management helps produce simplified, predictable, and beautiful results
- The polished collar on Straumann tissue level implants is designed to help soft tissue management
- Consistent Emergence Profiles™ in Straumann bone level implants help ensure tissue is managed the same way throughout the restoration: from healing cap to crown, producing the right esthetic result



Straumann CAD/CAM abutments

- Custom milled zirconia dioxide or titanium abutments are designed for outstanding esthetic results




Photos courtesy of Dr. Steven Fuqua, Southlake, Texas



Photos courtesy of Dr. Jason Gillespie, San Antonio, Texas

STRAUMANN® DENTAL IMPLANT SYSTEM

Experience confidence in implant treatment
with the Straumann Dental Implant System



- One system to meet the implant treatment needs of surgeon and referring dentist
- Reliable soft and hard tissue maintenance through multifactorial Bone Control Design™
- Faster osseointegration and increased predictability in stability-critical treatment protocols thanks to the unique SLActive® surface*
- Simple restorative options for a variety of implant treatments
- Peace of mind through high survival rates and extensive scientific evidence¹

*As compared to SLA® in animal studies

¹ Data on File



STRAUMANN GUARANTEE

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01

North American Distributors

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168
978/747 2500
Fax 978/747 2490
www.straumannusa.com

Straumann Canada Limited
4145 North Service Road, Suite 303
Burlington, ON L7L 6A3
Phone 800/363 4024
905/319 2900
Fax 905/319 2911
www.straumann.ca