

STRAUMANN SIMPLICITY

BONE LEVEL CEMENTABLE ABUTMENT ALL-IN-ONE SET

Straumann's bone level platform offers treatment flexibility and restorative simplicity for delivering esthetic results patients demand. A broad array of stock abutments and custom abutments offer clinicians a range of material, height, angulation and diameter options.

Ordering and stocking components is now easier than ever with Straumann's Bone Level Cementable Abutment All-In-One Sets. These convenient kits contain all the components you need for the restoration – cementable abutment, analog, protective cap, burnout coping, impression cap, and basal screw – in a single package for straightforward inventory and case planning.



Article Number	Article Description
027.4301	NC Cem. All-In-One Set, D3.5, AH4, GH1, crown
027.4302	NC Cem. All-In-One Set, D3.5, AH4, GH2, crown
027.4303	NC Cem. All-In-One Set, D3.5, AH4, GH3, crown
027.4307	NC Cem. All-In-One Set, D5, AH4, GH1, crown
027.4308	NC Cem. All-In-One Set, D5, AH4, GH2, crown
027.4309	NC Cem. All-In-One Set, D5, AH4, GH3, crown
027.4310	NC Cem. All-In-One Set, D5, AH5.5, GH1, crown
027.4311	NC Cem. All-In-One Set, D5, AH5.5, GH2, crown
027.4312	NC Cem. All-In-One Set, D5, AH5.5, GH3, crown

Article Number	Article Description
027.4313	RC Cem. All-In-One Set, D5, AH4, GH1, crown
027.4314	RC Cem. All-In-One Set, D5, AH4, GH2, crown
027.4315	RC Cem. All-In-One Set, D5, AH4, GH3, crown
027.4316	RC Cem. All-In-One Set, D5, AH5.5, GH1, crown
027.4317	RC Cem. All-In-One Set, D5, AH5.5, GH2, crown
027.4318	RC Cem. All-In-One Set, D5, AH5.5, GH3, crown
027.4319	RC Cem. All-In-One Set, D6.5, AH4, GH1, crown
027.4320	RC Cem. All-In-One Set, D6.5, AH4, GH2, crown
027.4321	RC Cem. All-In-One Set, D6.5, AH4, GH3, crown
027.4322	RC Cem. All-In-One Set, D6.5, AH5.5, GH1, crown
027.4323	RC Cem. All-In-One Set, D6.5, AH5.5, GH2, crown
027.4324	RC Cem. All-In-One Set, D6.5, AH5.5, GH3, crown

For additional information or to place an order, please contact Straumann Customer Service at 800/448 8168 or your local Straumann Territory Manager directly.