

More than a barrier.

Consistent support.

Straumann® Membrane Plus™ enhances graft stabilization and bone growth by **providing soft tissue support and space maintenance over a predictable timeframe.** Meticulously manufactured from bovine material that is one of the purest forms of collagen, it's treated and sterilized to be biocompatible, and crosslinked for predictable resorption rates. With its easy handling characteristics, it's designed to be the easy way to provide critical support.

EFFECTIVE REINFORCEMENT

- Cell-occlusive to support bone regeneration, protecting the graft area from unwanted soft tissue infiltration during the initial phase of healing
- Predictably resorbs after 6 months, eliminating the need for a second surgery for removal

DURABLE AND LONG LASTING


- Composed of multiple layers of type 1 collagen, providing soft tissue support
- Slow resorption over 6 months
- Crosslinking allows for reliable performance even if exposed

FLEXIBILITY TO MEET CLINICAL NEEDS


- Combine with Straumann AlloGraft and/or Straumann XenoGraft to maintain ideal space and long-term cell occlusion for maximum bone volume
- Provides support for peri-implant dehiscence defects, fenestration defects, sinus lifts, and horizontal or ridge augmentations

EASY TO HANDLE


- Rigid membrane, yet only 0.3 mm thick
- Drapability facilitates ease of use
- Not side-specific, enabling easy placement
- Minimal hydration creates optimal bio-adaptability
- Minimal suturing needed, potentially reducing procedure time and patient discomfort




PROCESSING ASSURES COMPATIBILITY



CASE STUDY



Site grafted with xenograft



Barrier membrane in place



Primary closure achieved



Three-month healing



Final prosthesis in place

Photos courtesy of Dr. Sungmin Chang

To learn more, **contact your Straumann Territory Manager or Customer Service.**

Buy online at straumann.us/eShop.



International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone: +41 (0)61 965 11 11
Fax: +41 (0)61 965 11 01

Straumann North American Headquarters

Straumann USA, LLC www.straumann.us
60 Minuteman Road www.straumann.ca
Andover, MA 01810
Phone: 800.448.8168 (US) • 800.363.4024 (CA)
Fax: 978.747.2490

© Straumann USA, LLC 2014. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved. May not be available in all countries, subject to regulatory approval and market release. Please ask your Straumann Territory Manager for more information and availability.