

More than Tissue Regeneration. The Natural Choice.

Straumann® AlloGraft Dermal Matrix is an acellular dermal matrix excellent for GTR/GBR and oral soft tissue correction in areas such as periodontal defects and ridge preservation.

BENEFITS OF STRAUMANN® DERMAL MATRIX:


- **Safety:** Packaged Sterile with an SAL of 10^{-6} for increased patient safety
- **Convenient Storage:** Ambient temperature, no refrigeration is needed
- **Biocompatible:** $\geq 97\%$ of DNA removed while retaining native growth factors, collagen and elastin
- **Ready to use:** No rehydration needed
- **Protects Graft Site:** Cell-occlusive natural membrane that protects grafted site from fluid and bacterial contamination
- **Effective Site Remodeling:** Graft resorbs through normal site remodeling in 4-6 months following placement

CLINICAL APPLICATIONS:

- Membrane (Bone Graft Containment)
- Gingival Recession
- Implant preparation
- Ridge Preservation
- Periodontal Defects
- Soft Tissue Repair/Correction



PROCESSING ASSURES COMPATIBILITY



Product Specifications	
Article Number	Product Thickness 0.76-1.25 mm
DM-100	15x20 mm
DM-101	20x40 mm
DM-150	10x10 mm
DM-151	10x40 mm

Product Specifications	
Article Number	Product Thickness 1.26-1.75 mm
DM-200	15x20 mm
DM-201	20x40 mm
DM-250	10x10 mm
DM-251	10x40 mm

To learn more, **contact your Straumann Territory Manager or Customer Service.**

Buy online at straumann.us/eShop or straumann.ca/eShop.



Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168 (US) • 800/363 4024 (CA)
Fax 978/747 2490
www.straumann.us • www.straumann.ca

© Straumann USA, LLC 2017. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.

LifeNet Health® and Matracell® are registered trademarks of LifeNet Health, Inc., Virginia Beach, VA.