



# LOCATOR®


## Restorative Components


### BONE LEVEL IMPLANTS


Article Number		Article
<b>NC Bone Level LOCATOR Abutments</b>		
Restorative platform Ø 3.3		
022.2502		NC LOCATOR, H 2.0 mm
022.2503		NC LOCATOR, H 3.0 mm
022.2504		NC LOCATOR, H 4.0 mm
022.2505		NC LOCATOR, H 5.0 mm
022.2506		NC LOCATOR, H 6.0 mm

<b>RC Bone Level LOCATOR Abutments</b>		
Restorative platform Ø 4.1 or Ø 4.8		
022.4501		RC LOCATOR, H 1.0 mm
022.4502		RC LOCATOR, H 2.0 mm
022.4503		RC LOCATOR, H 3.0 mm
022.4504		RC LOCATOR, H 4.0 mm
022.4505		RC LOCATOR, H 5.0 mm
022.4506		RC LOCATOR, H 6.0 mm

### TISSUE LEVEL IMPLANTS

Article Number		Article
<b>RN Tissue Level LOCATOR Abutments</b>		
Restorative platform Ø 4.8		
048.175		RN LOCATOR, H 1.0 mm
048.176		RN LOCATOR, H 2.0 mm
048.177		RN LOCATOR, H 3.0 mm
048.178		RN LOCATOR, H 4.0 mm
048.179		RN LOCATOR, H 5.0 mm
048.180		RN LOCATOR, H 6.0 mm

<b>WN Tissue Level LOCATOR Abutments</b>		
Restorative platform Ø 6.5		
048.183		WN LOCATOR, H 1.0 mm
048.184		WN LOCATOR, H 2.0 mm
048.185		WN LOCATOR, H 3.0 mm
048.186		WN LOCATOR, H 4.0 mm
048.187		WN LOCATOR, H 5.0 mm

<b>NNC Tissue Level LOCATOR Abutments</b>		
Restorative platform Ø 3.5		
045.581		NNC LOCATOR, H 1.0 mm
045.582		NNC LOCATOR, H 2.0 mm
045.583		NNC LOCATOR, H 3.0 mm
045.584		NNC LOCATOR, H 4.0 mm
045.585		NNC LOCATOR, H 5.0 mm
045.586		NNC LOCATOR, H 6.0 mm






**Additional components may be available.**




**For more information or to place an order, contact your local Straumann representative or Customer Service at 800/448 8168 (US) or 800/363 4024 (CA).**

**TIP:** To determine correct locator cuff height, we recommend measuring from the top of the implant to the highest point of tissue (rounding up).












H = Height  
 NNC = Narrow Neck Crossfit®  
 RN = Regular Neck  
 WN = Wide Neck  
 NC = Narrow Crossfit  
 RC = Regular Crossfit

## AUXILLARY TOOLS & PARTS

Article Number		Article
046.415		LOCATOR Core Tool
046.416		LOCATOR Driver, short for ratchet, L 15.0 mm
046.417		LOCATOR Driver, Long for ratchet, L 21.0 mm
048.199v4		LOCATOR parallel post
048.200		LOCATOR angle measurement guide
046.414		LOCATOR Core Tool tip

Article Number		Article
<b>Impression Components</b>		
025.0120-04		BL LOCATOR Female Analog
048.198v4		TL LOCATOR Female Analog
048.197v4		LOCATOR Impression Coping

"-04" "v4" = components come in packs of 4

Article Number		Article	Dimensions	Material
<b>Universal Male Processing Parts</b>				
048.189V2		LOCATOR male processing package, contents: denture cap (∅ 5.5 mm, height 2.5 mm) with black processing male (height 1.9 mm), block-out spacer (048.196), white replacement male (048.190), pink replacement male (048.191), blue replacement male (048.192)		titanium/nylon
048.190V4		LOCATOR replacement male, white, 0°–10°, 5 lbs, 2.27 kg**	height 1.7 mm	nylon
048.191V4		LOCATOR replacement male, pink, light retention, 0°–10°, 3 lbs, 1.36 kg**	height 1.7 mm	nylon
048.192V4		LOCATOR replacement male, blue, extra-light retention, 0°–10°, 1.5 lbs, 0.68 kg**	height 1.7 mm	nylon
048.182V2		LOCATOR male processing package, extended range, contents: denture cap (∅ 5.5 mm, height 2.5 mm) with black processing male (height 1.9 mm), block-out spacer (048.196), green replacement male (048.193), orange replacement male (048.188), red replacement male (048.194)		titanium/nylon
048.193V4		LOCATOR replacement male, green, extended range, 10°–20°, 3–4 lbs, 1.36–1.82 kg**	height 1.7 mm	nylon
048.188V4		LOCATOR replacement male, orange, light retention, extended range, 10°–20°, 2 lbs, 0.91 kg**	height 1.7 mm	nylon
048.194V4		LOCATOR replacement male, red, extra-light retention, extended range, 10°–20°, 0.5–1.5 lbs, 0.23–0.68 kg**	height 1.7 mm	nylon
048.181V4		LOCATOR replacement male, grey, extra-light retention, extended range, 0 lbs.	height 1.7 mm	nylon
048.218V4		LOCATOR processing spacer		POM
046.413V4		LOCATOR abutment holder sleeve		PSU

\* = For the correction of angle divergences

\*\* = Retention force