








Simplified healing with the Straumann® bone level implant.

BONE LEVEL IMPLANT DIAMETERS

Narrow CrossFit® Connection = NC*	Regular CrossFit Connection = RC*	
 Ø 3.3 mm	 Ø 4.1 mm	 Ø 4.8 mm

* NC and RC implants are available in lengths of 8.0 mm, 10.0 mm, 12.0 mm, and 14.0 mm

































CLOSURE SCREWS

Narrow CrossFit Connection		Regular CrossFit Connection	
 0.0 mm	 0.5 mm	 0.0 mm	 0.5 mm

HEALING ABUTMENTS

Diameter	Ø 3.6 mm	Ø 4.8 mm	Ø 3.3 mm	Ø 5.0 mm	Diameter	Ø 4.5 / 5.0 mm	Ø 6.0 / 6.5 mm	Ø 4.4 / 4.7 mm	Ø 7.0 mm
									
Available heights	5.0 mm 3.5 mm 2.0 mm	5.0 mm 3.5 mm 2.0 mm	5.0 mm 3.5 mm			6.0 mm 4.0 mm 2.0 mm	6.0 mm 4.0 mm 2.0 mm	6.0 mm 4.0 mm	

OVERVIEW OF CLOSURE SCREWS AND HEALING ABUTMENTS

Narrow CrossFit® Connection = NC				
Connection		Article Description		Article No.
Submucosal healing	NC		NC Closure cap, H 0	024.2100S
	NC		NC Closure cap, H 0 (pack of 4)	024.2100-04†
	NC		NC Closure cap, H 0.5	024.2105S
	NC		NC Closure cap, H 0.5 (pack of 4)	024.2105-04†
	NC		NC Healing Abt., conical, D 3.6, H 2.0	024.2222S
Transmucosal healing	NC		NC Healing Abt., conical, D 3.6, H 3.5	024.2224S
	NC		NC Healing Abt., conical, D 3.6, H 5.0	024.2226S
	NC		NC Healing Abt., conical, D 4.8, H 2.0	024.2242S
	NC		NC Healing Abt., conical, D 4.8, H 3.5	024.2244S
	NC		NC Healing Abt., conical, D 4.8, H 5.0	024.2246S
Esthetic region	NC		NC Healing Abt., bottle, D 3.3, H 3.5	024.2234S
	NC		NC Healing Abt., bottle, D 3.3, H 5.0	024.2236S
	NC		NC Healing Abt., custom, D 5.0	024.2270S
Regular CrossFit Connection = RC				
Connection		Article Description		Article No.
Submucosal healing	RC		RC Closure cap, H 0	024.4100S
	RC		RC Closure cap, H 0 (pack of 4)	024.4100-04†
	RC		RC Closure cap, H 0.5	024.4105S
	RC		RC Closure cap, H 0.5 (pack of 4)	024.4105-04†
	RC		RC Healing Abt., conical, D 4.5, H 2.0	024.0000S
	RC		RC Healing Abt., conical, D 4.5, H 4.0	024.0001S
	RC		RC Healing Abt., conical, D 4.5, H 6.0	024.0002S
	RC		RC Healing Abt., conical, D 6.0, H 2.0	024.0003S
	RC		RC Healing Abt., conical, D 6.0, H 4.0	024.0004S
	RC		RC Healing Abt., conical, D 6.0, H 6.0	024.0005S
	RC		RC Healing Abt., conical, D 5.0, H 2.0	024.4222S*
Transmucosal healing	RC		RC Healing Abt., conical, D 5.0, H 4.0	024.4224S*
	RC		RC Healing Abt., conical, D 5.0, H 6.0	024.4226S*
	RC		RC Healing Abt., conical, D 6.5, H 2.0	024.4242S*
	RC		RC Healing Abt., conical, D 6.5, H 4.0	024.4244S*
	RC		RC Healing Abt., conical, D 6.5, H 6.0	024.4246S*
Esthetic region	RC		RC Healing Abt., bottle, D 4.4, H 4.0	024.4234S
	RC		RC Healing Abt., bottle, D 4.7, H 6.0	024.4236S
	RC		RC Healing Abt., custom, D 7.0	024.4270S

S = Sterile

H = Height

D = Diameter

† Sterile pack of 4

* Indicates conical BL RC Healing Abutment with improved shape

All measurements in millimeters unless otherwise specified.