



ESTHETICS - BY DESIGN

**STRAUMANN® ANATOMIC
IPS E.MAX® ABUTMENT**
BY STRAUMANN
AND IVOCLAR VIVADENT®



COMMITTED TO
SIMPLY DOING MORE
FOR DENTAL PROFESSIONALSSM

ESTHETICS -

STRAUMANN® ANATOMIC IPS E.MAX® ABUTMENT

Straumann is the partner of choice for innovations in implant and restorative dentistry. Ivoclar Vivadent is an established specialist in ceramic materials and final restorations. The synergy of technologies has resulted in a premium esthetic solution.



STRAUMANN AND IVOCALAR VIVADENT PRESENT THE STRAUMANN ANATOMIC IPS E.MAX® ABUTMENT

- Designed for both the NC and RC Bone Level Connection
- Highly esthetic ZrO₂ abutment in multiple shades for natural appearance
- Designed for easy-to-achieve esthetics and strength
- Designed and precision made by Straumann and Ivoclar Vivadent

FOR ALL-CERAMIC CROWNS FROM ABUTMENT TO FINAL RESTORATION



BY DESIGN

ENT BY STRAUMANN AND IVOCCLAR VIVADENT®



HIGH-END ESTHETICS

- Pre-shaped abutments with prepared mucosa margins helps in creating efficient restorations
- Eliminates visible metal margins, creating natural smiles
- Highly esthetic results possible even with thin tissue biotypes

FLEXIBILITY

- Designed to modify it chairside
- A variety of treatment options
- Available in multiple shades for patient needs

PREDICTABILITY

- Deliver reliable, natural looking results
- Precision-made pre-fabricated abutments for consistent results

Color





straight

15° angled

MO 0 (white)

NC Article Number





RC Article Number

			
022.2812	022.2814	022.2822	022.2824
022.4812	022.4814	022.4822	022.4824

MO1 (shaded)

NC Article Number

RC Article Number

			
022.2832	022.2834	022.2842	022.2844
022.4832	022.4834	022.4842	022.4844

GH 2mm

GH 3.5mm

GH 2mm

GH 3.5mm



COMBINED EXPERTISE - VISIBLE RESULTS

HIGH-END ESTHETICS – all-ceramic material for natural appearance

FLEXIBILITY – a multitude of restorative treatment options

PREDICTABILITY – complementary features designed for beautiful results

**CONTACT STRAUMANN CUSTOMER SERVICE AT 800/448 8168
TO LEARN MORE ABOUT STRAUMANN ESTHETICS OR CONTACT
A REPRESENTATIVE IN YOUR AREA.**

www.straumann.com