




## ACHIEVE MORE with **coDiagnostiX™ 8.0**



**3D**



### IMPROVE YOUR PRACTICE BRAND WITH NEW MARKETING TOOLS

- Present to your patients and referrals with the free iPad® application from the Apple® iPad App Store
- Improve patient case acceptance by sharing treatment plan visualization in real 3D
- Request expert opinions through the free coDiagnostiX caseXchange™ server

### COLLABORATE

- Share your treatment plans via iPad with coDiagnostiX Presenter
- Consult with other coDiagnostiX users with CARES® caseXchange™
- Visualize your treatment plan in 3 dimensions


### EXPAND YOUR CAPABILITIES

- Import STL files from iTero® or other scanners
- Measure distances between cemento-enamel junction and the bone level with coPeriodontiX™
- Design your restorative plan with a Virtual Tooth Library
- Align dataset conveniently in 3 dimensions
- Place implants directly through the surgical guide with Straumann® Guided Implants

## A POWERFUL TOOL FOR PROSTHETIC DRIVEN **IMPLANT PLANNING**



COMMITTED TO  
**SIMPLY DOING MORE**  
FOR DENTAL PROFESSIONALS<sup>SM</sup>



**NEW STRAUMANN GUIDED IMPLANTS - FOR FULLY GUIDED SURGERY\***

- Fully guided implant insertion through surgical templates containing Straumann Ø 5.0 mm sleeves
- Guided implant mount is pre-mounted and is firmly attached to the implant
- Physical depth control with stop key
- Visual depth control
- Help to perform minimally invasive surgery with new mucosa punches

\* Please refer to the Straumann Product Catalog (USLIT 007) for complete portfolio offerings for Straumann Guided Implants.



**NEW Pay Per Guide software option** for unlimited treatment planning  
Additional fee for surgical template planning

**PC Visit remote support** allows your Straumann technical expert to log on to your computer and collaborate in real time

**To learn more about Straumann Guided Surgery, please contact your local sales representative or Straumann Customer Service at 800/448 8168.**

[www.straumann.com](http://www.straumann.com)



STRAUMANN GUARANTEE

**International Headquarters**

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
Fax +41 (0)61 965 11 01

**Straumann USA**

Straumann USA, LLC  
60 Minuteman Road  
Andover, MA 01810  
Phone 800/448 8168  
978/747 2500  
Fax 978/747 2490  
[www.straumannusa.com](http://www.straumannusa.com)

**Straumann Canada**

Straumann Canada Limited  
3115 Harvester Road, Suite 100  
Burlington, ON L7N 3N8  
Phone 800/363 4024  
905/319 2900  
Fax 905/319 2911  
[www.straumann.ca](http://www.straumann.ca)

© Straumann USA, LLC 2011. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved. iPad® is a trademark of Apple Inc., registered in the U.S. and other countries PC Visit is a product of PC Visit Software AG.