

## IMPLANTS FOR THE STRAUMANN® DENTAL IMPLANT SYSTEM

Customer #: \_\_\_\_\_ Name: \_\_\_\_\_ Date: \_\_\_\_\_


Payment type: NET 30 Credit Card Fed Ex: Overnight/10:30am Overnight/4:30pm 2-Day 3-5 Day Page \_\_\_\_\_ of \_\_\_\_\_

Enter quantities in first column of each row. Send completed fax form to 800/524 6752.

Standard Ø 3.3 mm, RN		
033.431S		Roxolid® 8.0 mm
033.131S		SLActive 8.0 mm
043.131S		SLA 8.0 mm
033.432S		Roxolid 10.0 mm
033.132S		SLActive 10.0 mm
043.132S		SLA 10.0 mm
033.433S		Roxolid 12.0 mm
033.133S		SLActive 12.0 mm
043.133S		SLA 12.0 mm
033.434S		Roxolid 14.0 mm
033.134S		SLActive 14.0 mm
043.134S		SLA 14.0 mm
033.435S		Roxolid 16.0 mm
033.135S		SLActive 16.0 mm
043.135S		SLA 16.0 mm


Standard Ø 4.1 mm, RN		
033.030S		SLActive 6.0 mm
043.030S		SLA 6.0 mm
033.031S		SLActive 8.0 mm
043.031S		SLA 8.0 mm
033.032S		SLActive 10.0 mm
043.032S		SLA 10.0 mm
033.033S		SLActive 12.0 mm
043.033S		SLA 12.0 mm
033.034S		SLActive 14.0 mm
043.034S		SLA 14.0 mm
033.035S		SLActive 16.0 mm
043.035S		SLA 16.0 mm


Standard Ø 4.8 mm, RN		
033.230S		SLActive 6.0 mm
043.230S		SLA 6.0 mm
033.231S		SLActive 8.0 mm
043.231S		SLA 8.0 mm
033.232S		SLActive 10.0 mm
043.232S		SLA 10.0 mm
033.233S		SLActive 12.0 mm
043.233S		SLA 12.0 mm
033.234S		SLActive 14.0 mm
043.234S		SLA 14.0 mm


Standard Ø 4.8 mm, WN		
033.630S		SLActive 6.0 mm
043.630S		SLA 6.0 mm
033.631S		SLActive 8.0 mm
043.631S		SLA 8.0 mm
033.632S		SLActive 10.0 mm
043.632S		SLA 10.0 mm
033.633S		SLActive 12.0 mm
043.633S		SLA 12.0 mm

Standard Plus Ø 3.3 mm, NN		
033.951S		SLActive 8.0 mm
042.930S		SLA 8.0 mm
033.952S		SLActive 10.0 mm
042.931S		SLA 10.0 mm
033.953S		SLActive 12.0 mm
042.932S		SLA 12.0 mm
033.954S		SLActive 14.0 mm
042.933S		SLA 14.0 mm


Standard Plus Ø 3.3 mm, RN		
033.451S		Roxolid 8.0 mm
033.151S		SLActive 8.0 mm
043.151S		SLA 8.0 mm
033.452S		Roxolid 10.0 mm
033.152S		SLActive 10.0 mm
043.152S		SLA 10.0 mm
033.453S		Roxolid 12.0 mm
033.153S		SLActive 12.0 mm
043.153S		SLA 12.0 mm
033.454S		Roxolid 14.0 mm
033.154S		SLActive 14.0 mm
043.154S		SLA 14.0 mm


Standard Plus Ø 4.1 mm, RN		
033.050S		SLActive 6.0 mm
043.050S		SLA 6.0 mm
033.051S		SLActive 8.0 mm
043.051S		SLA 8.0 mm
033.052S		SLActive 10.0 mm
043.052S		SLA 10.0 mm
033.053S		SLActive 12.0 mm
043.053S		SLA 12.0 mm
033.054S		SLActive 14.0 mm
043.054S		SLA 14.0 mm

Standard Plus Ø 4.8 mm, RN		
033.250S		SLActive 6.0 mm
043.250S		SLA 6.0 mm
033.251S		SLActive 8.0 mm
043.251S		SLA 8.0 mm
033.252S		SLActive 10.0 mm
043.252S		SLA 10.0 mm
033.253S		SLActive 12.0 mm
043.253S		SLA 12.0 mm
033.254S		SLActive 14.0 mm
043.254S		SLA 14.0 mm

Standard Plus Ø 4.8 mm, WN		
033.650S		SLActive 6.0 mm
043.650S		SLA 6.0 mm
033.651S		SLActive 8.0 mm
043.651S		SLA 8.0 mm
033.652S		SLActive 10.0 mm
043.652S		SLA 10.0 mm
033.653S		SLActive 12.0 mm
043.653S		SLA 12.0 mm

Tapered Effect Ø 3.3 mm, RN		
033.751S		Roxolid 8.0 mm
033.721S		SLActive 8.0 mm
043.721S		SLA 8.0 mm
033.752S		Roxolid 10.0 mm
033.722S		SLActive 10.0 mm
043.722S		SLA 10.0 mm
033.753S		Roxolid 12.0 mm
033.723S		SLActive 12.0 mm
043.723S		SLA 12.0 mm
033.754S		Roxolid 14.0 mm
033.724S		SLActive 14.0 mm
043.724S		SLA 14.0 mm

Tapered Effect Ø 4.1 mm, RN		
033.761S		SLActive 8.0 mm
043.761S		SLA 8.0 mm
033.762S		SLActive 10.0 mm
043.762S		SLA 10.0 mm
033.763S		SLActive 12.0 mm
043.763S		SLA 12.0 mm
033.764S		SLActive 14.0 mm
043.764S		SLA 14.0 mm

Tapered Effect Ø 4.8 mm, WN		
033.712S		SLActive 10.0 mm
043.712S		SLA 10.0 mm
033.713S		SLActive 12.0 mm
043.713S		SLA 12.0 mm
033.714S		SLActive 14.0 mm
043.714S		SLA 14.0 mm

**Important:**

Be sure to select the correct implant surface (SLActive® or SLA®) when ordering. All Roxolid implants feature the SLActive® surface.




## IMPLANTS FOR THE STRAUMANN® DENTAL IMPLANT SYSTEM
















Customer #: \_\_\_\_\_ Name: \_\_\_\_\_ Date: \_\_\_\_\_

Payment type: NET 30 Credit Card Fed Ex: Overnight/10:30am Overnight/4:30pm 2-Day 3-5 Day Page \_\_\_\_\_ of \_\_\_\_\_



Enter quantities in first column of each row. Send completed fax form to 800/524 6752.

### Tissue Level Healing Caps

Straumann® Bone Level Implant, Ø 3.3 mm, NC		
021.2208		Roxolid® 8.0 mm
021.2108		SLActive 8.0 mm
021.2408		SLA 8.0 mm
021.2210		Roxolid® 10.0 mm
021.2110		SLActive 10.0 mm
021.2410		SLA 10.0 mm
021.2212		Roxolid® 12.0 mm
021.2112		SLActive 12.0 mm
021.2412		SLA 12.0 mm
021.2214		Roxolid® 14.0 mm
021.2114		SLActive 14.0 mm
021.2414		SLA 14.0 mm
Ø 4.1 mm, RC		
021.4108		SLActive 8.0 mm
021.4408		SLA 8.0 mm
021.4110		SLActive 10.0 mm
021.4410		SLA 10.0 mm
021.4112		SLActive 12.0 mm
021.4412		SLA 12.0 mm
021.4114		SLActive 14.0 mm
021.4414		SLA 14.0 mm
Ø 4.8 mm, RC		
021.6108		SLActive 8.0 mm
021.6408		SLA 8.0 mm
021.6110		SLActive 10.0 mm
021.6410		SLA 10.0 mm
021.6112		SLActive 12.0 mm
021.6412		SLA 12.0 mm
021.6114		SLActive 14.0 mm
021.6414		SLA 14.0 mm

NN healing caps		
048.374		H 1.8 mm, Ti
048.043		H 3.4 mm, Ti
048.050		3.4 mm, PEEK
RN closure screws and healing caps		
048.371V4		Internal
048.373V4		H 1.5 mm
048.033		H 2.0 mm
048.034		H 3.0 mm
048.037		H 4.5 mm
048.028		H 2.0 mm, beveled
048.029		H 3.5 mm, beveled
WN closure screw and healing caps		
048.375		Internal
048.038		H 2.0 mm
048.039		H 3.0 mm
048.053		H 4.5 mm
048.030		H 2.0 mm, beveled

### Bone Level Healing Caps

Healing Abutment, customizable*		
024.2270		NC, D 5.0 mm
024.4270		RC, D 7.0 mm
Healing Abutment, conical shape		
024.2222		NC, D 3.6 mm, GH 2.0 mm
024.2224		NC, D 3.6 mm, GH 3.5 mm
024.2226		NC, D 3.6 mm, GH 5.0 mm
024.2242		NC, D 4.8 mm, GH 2.0 mm
024.2244		NC, D 4.8 mm, GH 3.5 mm
024.2246		NC, D 4.8 mm, GH 5.0 mm
024.4222		RC, D 4.5 mm, GH 2.0 mm
024.4224		RC, D 4.5 mm, GH 4.0 mm
024.4226		RC, D 4.5 mm, GH 6.0 mm
024.4242		RC, D 6.0 mm, GH 2.0 mm
024.4244		RC, D 6.0 mm, GH 4.0 mm
024.4246		RC, D 6.0 mm, GH 6.0 mm
Healing Abutment, bottle shape		
024.2234		NC, D 3.3 mm, GH 3.5 mm
024.2236		NC, D 3.3 mm, GH 5.0 mm
024.4234		RC, D 4.4 mm, GH 4.0 mm
024.4236		RC, D 4.7 mm, GH 6.0 mm

V4 = 4 components per pack. All measurements given are in millimeters unless otherwise specified.

- SLA®** Sand-blasted, large grit, Acid-etched
  - SLActive®** This roughened implant surface shares the same characteristics and topography as SLA® implants with the additional properties of chemical activity and hydrophilicity
  - Roxolid®** Titanium Zirconium alloy implant material
  - S** Standard (implants that have a 2.8 mm smooth collar height)
  - SP** Standard Plus (implants that have a 1.8 mm smooth collar height)
  - TE** Tapered Effect (implants that combine a cylindrical-shaped apical portion with a tapered section in the conical endosteal region. They also have a 1.8 mm smooth collar height)
  - NN** Narrow Neck (Ø 3.5 mm restorative platform)
  - RN** Regular Neck (Ø 4.8 mm restorative platform)
  - WN** Wide Neck (Ø 6.5 mm restorative platform)
  - NC** Narrow CrossFit® Connection
  - RC** Regular CrossFit® Connection
- \* corresponding basal screw included in package