

The next generation in drill technology



Straumann® Twist Drill PRO Ø 4.2 mm

Careful surgical site planning is critical for obtaining high implant primary stability. For optimal implant bed preparation, drills with wide diameters require particularly precise handling.

The Straumann® Twist Drill PRO Ø 4.2 mm is designed to help prevent over-preparation of the osteotomy site. The rounded flat-end tip facilitates safe drill positioning and ensures secure drill centering in the implant bed. The defined direction is maintained throughout the drilling process. The Straumann® Twist Drill PRO Ø 4.2 mm is available in lengths of 3.3 mm (short) and 4.1 mm (long).